

User's Guide




**THE CENTER FOR
NANOSCIENCE AND
NANOTECHNOLOGY**

MA6 MASK ALIGNER






Machine Overview





Part #	Description
1	Log book
2	<i>MAIN POWER SWITCH</i>
3	Lamp power supply
4	Analog panel
5	Digital panel
6	Microscope
7	Mask holder

Part #	Description
8	Mask holder vacuum connector
9	Slide
10	Chuck
11	Chuck positioning micrometers θ (front), x-axis (rear)
12	Chuck positioning micrometer y-axis
13	Shutter




Lamp Start and Warm Up

1	Log in on Log Book	
2	<p style="text-align: center;">Turn <i>MAIN POWER SWITCH</i> CLOCKWISE</p> <p><i>ON</i> displayed in small window a</p> <p>Power supply will read <i>STAND-BY</i></p>	 <p style="text-align: center;"><i>Part # 1</i></p>
3	<p>a Wait about 10 seconds and press <i>ON</i> on the power supply</p> <p>Power supply will read </p> <p>b Wait about 10 seconds and press <i>CP</i> on the power supply</p> <p>Power supply will read </p> <p>c Wait about 10 seconds and press <i>START</i> on the power supply</p> <p>Power supply will read </p>	 <p style="text-align: center;"><i>Part # 2</i></p>
4	<p>After several minutes values will appear on screen</p> <p>Wait about 15 minutes before turning machine on</p>	
*	<p>If power supply beeps and displays <i>VOLT ERROR</i></p> <p style="text-align: center;">Press <i>OFF</i> and then <i>ON</i> on the power supply</p> <p>Wait until screen displays <i>STAND BY</i></p> <p>Perform STEP 3 again</p>	



Machine Start

<p>5</p>	<p>Turn machine on</p> <p>Turn machine power switch CLOCKWISE</p> <p>Switch will spring back to middle</p>	
<p>6</p>	<p>Wait until screen reads <i>Ready for Start</i></p> <p>Press LOAD button</p> <p>Wait until main screen is displayed</p>	
<p>7</p>	<p>a Contact <i>PROGRAM</i> Vacuum, Low-Vacuum, Hard, Soft and Proximity</p> <p>b Alignment type BSA—backside alignment (default) TSA—topside alignment</p> <p>c Alignment gap (in microns)</p> <p>d Microscope coordinates</p>	<p>Main Screen</p> 
<p>8</p>	<p>Replace chuck (if necessary)</p> <ol style="list-style-type: none"> 1. Gently pull slide out of machine until stop 2. Raise chuck from underneath 3. Switch chuck with one from shelf 4. Lower new chuck from underneath 5. Align white stripe with pin a 6. Gently push slide into machine 7. Align white crosses b with θ micrometer c 5. Realign a if necessary 	




Mask Insertion I

<p>9</p>	<p>Replacing mask holder (if necessary)</p> <p>Push red release ring towards machine</p> <p>Grasp metal connector on vacuum tube and extract from machine</p>	
<p>10</p>	<p>Do NOT touch the spacer flags! a</p> <p>Switch mask holder on side tray with one on shelf, holding at b and c</p> <p>Place new mask holder on side tray with spacer flags facing UP</p> <p>Grasp metal connector on vacuum tube and insert into machine</p>	
<p>11</p>	<p>Press <i>CHANGE MASK</i></p> <p>Mask vacuum is <i>OFF</i> by default</p>	
<p>12</p>	<ol style="list-style-type: none"> Place mask on holder, chrome side UP Align mask to left pin and top bar Push down metal plate to engage holder clip a 	






Mask Insertion II

<p style="font-size: 2em; font-weight: bold;">13</p>	<p style="text-align: center;">Press <i>ENTER</i> To engage mask vacuum</p> <p style="text-align: center;">Screen should read <i>ON</i></p> <p style="color: red;">Do not flip mask holder if mask is not held securely by vacuum</p>	 <p>The image shows a close-up of a control panel with a green 'ENTER' button. Below it is a green LCD screen displaying the following text: 'Hard-Ct/Cont. X :+ 0.00', 'TSA/A1. Gap: 45 Y :- 0.03', '- Change Mask - Press ENTER to toggle', and 'mask vacuum : ON'. A blue arrow points to the 'ENTER' button on the screen.</p>
<p style="font-size: 2em; font-weight: bold;">14</p>	<p style="text-align: center;">Grasp mask holder at b and c</p> <p style="text-align: center;">Flip mask holder over</p>	 <p>The image shows a gloved hand flipping a mask holder. Two blue circles labeled 'b' and 'c' indicate the points of contact for grasping the holder.</p>
<p style="font-size: 2em; font-weight: bold;">15</p>	<p style="text-align: center;">Align mask holder edge with rails a</p> <p style="color: red;">Gently slide mask holder into machine until stop</p>	 <p>The top image shows a gloved hand aligning the mask holder with a rail, indicated by a blue circle 'a' and an arrow. The bottom image shows the mask holder fully inserted into the machine.</p>
<p style="font-size: 2em; font-weight: bold;">16</p>	<p style="text-align: center;">Press <i>CHANGE MASK</i></p> <p style="text-align: center;">Engages vacuum between mask holder and machine</p> <p style="text-align: center;">Exits out of <i>CHANGE MASK</i> menu</p> <p style="text-align: center;">Returns to Main Screen</p>	 <p>The image shows a control panel with several buttons. A blue arrow points to the 'CHANGE MASK' button, which has a green indicator light above it. Other buttons include 'F1', 'LAMP TEST', 'SELECT PROGRAM', 'EDIT PROGRAM', 'EDIT PARAMETER', 'SCAN', 'SET REFERENCE', and 'ENTER'.</p>



Exposure without Alignment I

17	Main screen	
18	<p><u>To change Contact <i>PROGRAM</i></u></p> <p>Press <i>SELECT PROGRAM</i></p> <p>Use Y arrows on digital keypad to select program</p> <p>Press <i>SELECT PROGRAM</i> again to exit the menu</p>	 
19	<p><u>To change Exposure Time</u></p> <p>Press <i>EDIT PARAMETER</i></p> <p>Use Y arrows on digital keypad to change Exposure Time</p> <p>Press <i>EDIT PARAMETER</i> again to exit menu</p>	 
20	<p>Press <i>LOAD</i></p> <p>Gently pull slide out of machine until stop</p>	 

Exposure without Alignment II

<p style="font-size: 2em; font-weight: bold;">21</p>	<p>Place substrate on chuck and align to pins on 3", 4" and 6" chuck (4" chuck shown)</p> <p>or grooves on chip chuck (not shown)</p> <p style="color: red;">Do NOT let wafer sit on top of pin!</p> <p>When in doubt, leave a small gap (1 mm or so) between wafer and pins</p>	 <p>Shown: 4" wafer properly aligned to 1/2 pins</p> 
<p style="font-size: 2em; font-weight: bold;">22</p>	<p>When wafer is seated properly, press <i>ENTER</i></p> <p>To turn on substrate vacuum</p> <p style="color: red;">Gently push slide into machine until stop</p> <p>Press <i>ENTER</i> again</p> <p>Chuck raises automatically</p> <p>Screen reads <i>PERFORMING WEC</i></p>	  
<p style="font-size: 2em; font-weight: bold;">23</p>	<p style="color: red;">Make sure substrate is aligned to mask!</p> <p>If possible, align by looking through mask at substrate edges</p> <p>Otherwise, set X micrometer a and Y micrometer b both to 10 mm (center position)</p>	  
<p style="font-size: 2em; font-weight: bold;">24</p>	<p>Double check contact and exposure parameters!</p> <p style="color: red;">Put on UV protective glasses</p> <p>Press <i>EXPOSURE</i></p>	





Exposure without Alignment III

25	<p>Chuck will lower automatically</p> <p style="text-align: center;">When screen reads <i>Pull slide and unload expose substrate</i></p> <p style="text-align: center;">Gently pull slide until stop</p> <p>An audible beep is sounded when substrate vacuum is released</p>	
26	<p>Remove wafer from chuck</p>	
27	<p>Gently push slide into machine until stop</p>	
28	<p style="text-align: center;">Removing mask</p> <ol style="list-style-type: none"> 1. Press CHANGE MASK To release vacuum between mask holder and machine 2. Slide mask holder out of machine, flip and place on side tray spacer flags facing UP 3. Press ENTER To release substrate vacuum 5. Insert fingernail into groove a and pull TOWARDS YOU 6. Remove mask from mask holder 7. Press CHANGE MASK 8. Press ENTER To confirm mask holder is not in place, returns to main screen <p style="text-align: center;">During the day, leave MA6 in this state with Main Screen visible</p>	<div style="text-align: center;">  </div> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div>
29	<p>Log out on Log Book</p>	


Exposure with Alignment I

30	Main screen	
31	<p><u>To change Contact <i>PROGRAM</i></u></p> <p>Press <i>SELECT PROGRAM</i></p> <p>Use <i>Y</i> arrows on digital keypad to select program</p> <p>Press <i>SELECT PROGRAM</i> again to exit the menu</p>	
32	<p><u>To change Exposure Time</u></p> <p>Press <i>EDIT PARAMETER</i></p> <p>Use <i>Y</i> arrows on digital keypad to change Exposure Time</p> <p>Press <i>EDIT PARAMETER</i> again to exit menu</p>	
33	<p>Press <i>LOAD</i></p> <p>Gently pull slide out of machine until stop</p>	




Exposure with Alignment II

<p style="font-size: 2em; font-weight: bold;">34</p>	<p>Place substrate on chuck and align to pins on 3", 4" and 6" chuck (4" chuck shown)</p> <p>or grooves on chip chuck (not shown)</p> <p>Do NOT let wafer sit on top of pin!</p> <p>When in doubt, leave a small gap (1 mm or so) between wafer and pins</p>	 <p>Shown: 4" wafer properly aligned to 1/2 pins</p> 
<p style="font-size: 2em; font-weight: bold;">35</p>	<p>When wafer is seated properly,</p> <p style="text-align: center;">Press <i>ENTER</i></p> <p>To turn on substrate vacuum</p> <p>Gently push slide into machine until stop</p> <p style="text-align: center;">Press <i>ENTER</i> again</p> <p>Chuck raises automatically</p> <p>Screen reads <i>PERFORMING WEC</i></p>	  
<p style="font-size: 2em; font-weight: bold;">36</p>	<p style="text-align: center;">Press <i>BSA MICROSCOPE</i></p> <p>Turns off bottom side microscope</p> <p>Turns on top side microscope</p> <p><i>TSA displayed in main screen</i></p>	 
<p style="font-size: 2em; font-weight: bold;">37</p>	<p>Roughly align substrate to mask</p> <p>View substrate features through mask openings and use X (a) and Y (b) micrometers to align</p> <p>Otherwise set a and b to 10 mm</p>	 


Exposure with Alignment III

<p>38</p>	<p>To lower microscope</p> <p>Press <i>F1</i> and then <i>ENTER</i></p> <p>Caution!</p> <p>Make sure head is not under When it is lowering!</p>	 
<p>39</p>	<p>Turn <i>ILLUMINATION</i> knob to <i>TSA</i></p>	
<p>40</p>	<p>Rotate microscope objectives to lowest magnification (5X)</p>	
<p>41</p>	<p>Focus objectives with <i>TOP SUBSTRATE</i> fine focus knobs</p> <p>If focus is not possible, rotate microscope <i>Z</i> adjustment knob to the right of microscope and refocus with fine focus knobs</p>	
<p>42</p>	<p>Find alignment marks (AMs)</p> <p>Press <i>FAST</i> button a</p> <p>Use digital keypad b to move microscope and find AMs on one side</p> <p>Once found, use manual <i>X</i> adjusts c to find AMs on other side</p> <p>Use <i>SPLITFIELD</i> selector d as needed</p>	

Exposure with Alignment IV

43	<p>Check wafer/mask separation!</p> <p>Press <i>SEP</i> down button to view separation on screen</p> <p>While adjusting micrometers, minimum Z separation is 20 μm</p>	 
44	<p>Adjust X, Y and θ (a, b and c) micrometers to align substrate to mask</p>	
45	<p>Rotate microscope objectives to 10X to evaluate alignment (if necessary)</p>	
46	<p>Press <i>SEP</i> up button to raise chuck and evaluate alignment (if necessary)</p> <p>Caution!</p> <p>Do <u>NOT</u> adjust micrometers if Z separation is below 20 μm!</p> <p>If re-alignment is required Press <i>SEP</i> down button until screen reads 20 μm</p>	
47	<p>When alignment is achieved Double check contact and exposure parameters!</p> <p>Put on UV protective glasses</p> <p>Press <i>EXPOSURE</i></p>	 

Exposure with Alignment V

48	<p>Chuck will lower automatically</p> <p>When screen reads <i>Pull slide and unload expose substrate</i></p> <p style="color: red;">Gently pull slide until stop</p> <p>An audible beep is sounded when substrate vacuum is released</p>	
49	Remove wafer from chuck	
50	Gently push slide into machine until stop	
51	<p style="text-align: center;">Removing mask</p> <ol style="list-style-type: none"> 1. Press <i>CHANGE MASK</i> To release vacuum between mask holder and machine 2. Slide mask holder out of machine, flip and place on side tray spacer flags facing UP 3. Press <i>ENTER</i> To release substrate vacuum 5. Insert fingernail into groove and pull TOWARDS YOU 6. Remove mask from mask holder 7. Press <i>CHANGE MASK</i> 8. Press <i>ENTER</i> To confirm mask holder is not in place, returns to main screen <p style="color: red;">During the day, leave MA6 in this state with Main Screen visible</p>	  
52	Log out on Log Book	